

FIZ1 Antibody



PACO13438

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

FIZ1

Uniprot

Q96SL8

Synonyms:

FLT3-interacting zinc finger 1; FIZ1; FLJ00416; FLJ14768; ZNF798 ;

Immunogen:

Human domain-II-of-FIZ1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A