

PACO13449

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

EXOG

Uniprot

Q9Y2C4

Synonyms:

endo/exonuclease (5-3), endonuclease G-like; EXOG; ENDOGL1; ENDOGL2; ENGL; ENGL-B; ENGL-a; MGC125944; MGC125945 ;

Immunogen:

Human EXOG.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A