

MCM3AP Antibody



PACO13459

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

MCM3AP

Uniprot

O60318

Synonyms:

minichromosome maintenance complex component 3 associated protein; MCM3AP; FLJ44336; FLJ45306; GANP; KIAA0572; MAP80 ;

Immunogen:

Human GANP.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A