

GPRASP1 Antibody



PACO13468

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

GPRASP1

Uniprot

Q5JY77

Synonyms:

G protein-coupled receptor associated sorting protein 1; GPRASP1; GASP; GASP1; KIAA0443 ;

Immunogen:

Human GPRASP1-Specific.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A