

# GPRC5A Antibody



PACO13469

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

GPRC5A

**Uniprot**

Q8NFJ5

**Synonyms:**

G protein-coupled receptor, family C, group 5, member A; GPRC5A; GPCR5A; RAI3; RAIG1 ;

**Immunogen:**

Human GPRC5A, RAI3.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A