

GPRIN3 Antibody



PACO13470

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:**Protein Background:****Gene ID:**

GPRIN3

Uniprot

Q6ZVF9

Synonyms:

GPRIN family member 3; GPRIN3; FLJ42625; GRIN3 ;

Immunogen:

Human GPRIN3-Specific.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A