

BET1L Antibody



PACO13472

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

BET1L

Uniprot

Q9NYM9

Synonyms:

blocked early in transport 1 homolog (S. cerevisiae)-like; BET1L; BET1L1; GOLIM3; GS15; HSPC197 ;

Immunogen:

Human GS15.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A