

# GOSR1 Antibody



PACO13473

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

GOSR1

**Uniprot**

O95249

**Synonyms:**

golgi SNAP receptor complex member 1; GOSR1; GOLIM2; GOS-28; GOS28; GOS28/P28; GS28; P28 ;

**Immunogen:**

Human GS28.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A