

# ADH5 Antibody



PACO13474

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

ADH5

**Uniprot**

P11766

**Synonyms:**

alcohol dehydrogenase 5 (class III), chi polypeptide; ADH5; ADH-3; ADHX; FDH; GSNOR ;

**Immunogen:**

Human GSNOR, ADH5.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A