

HIST1H4A Antibody



PACO13479

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

HIST1H4A

Uniprot

P62805

Synonyms:

histone cluster 1, H4e; HIST1H4E; H4/j; H4FJ ;

Immunogen:

Human HIST1H4E.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A