

# LIN28B Antibody



PACO13509

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

LIN28B

**Uniprot**

Q6ZN17

**Synonyms:**

lin-28 homolog B (C. elegans); LIN28B; CSDD2; FLJ16517 ;

**Immunogen:**

Human Lin28B-specific.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A