

# LPAR3 Antibody



PACO13510

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

LPAR3

**Uniprot**

Q9UBY5

**Synonyms:**

lysophosphatidic acid, receptor 3; LPAR3; EDG7; Edg-7; FLJ98231; GPCR; HOFNH30; LPA3; LPA3; RP4-678I3 ;

**Immunogen:**

Human LPAR3-Specific.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A