

STARD3 Antibody



PACO13520

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

STARD3

Uniprot

Q14849

Synonyms:

StAR-related lipid transfer (START) domain containing 3; STARD3; CAB1; FLJ41370; MLN64; es64 ;

Immunogen:

Human MLN64, STARD3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A