

N6AMT1 Antibody



PACO13528

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

N6AMT1

Uniprot

Q9Y5N5

Synonyms:

N-6 adenine-specific DNA methyltransferase 1 (putative); N6AMT1; C21orf127; HEMK2; MGC19995; MTQ2; N6AMT; PRED28 ;

Immunogen:

Human N6AMT1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A