

SSNA1 Antibody



PACO13529

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SSNA1

Uniprot

O43805

Synonyms:

Sjogren syndrome nuclear autoantigen 1; SSNA1; N14; NA-14; NA14 ;

Immunogen:

Human NA14.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A