

NEO1 Antibody



PACO13537

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NEO1

Uniprot

Q92859

Synonyms:

neogenin homolog 1 (chicken); NEO1; DKFZp547A066; DKFZp547B146; HsT17534; IGDC2; NGN ;

Immunogen:

Human NEO1-Specific.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A