

SLCO1C1 Antibody



PACO13546

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SLCO1C1

Uniprot

Q9NYB5

Synonyms:

solute carrier organic anion transporter family, member 1C1; SLCO1C1; OATP-F; OATP1; OATP1C1; OATPRP5; SLC21A14 ;

Immunogen:

Human OATP14.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A