

RBL2 Antibody



PACO13550

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

RBL2

Uniprot

Q08999

Synonyms:

nucleolar and coiled-body phosphoprotein 1; NOLC1; KIAA0035; NOPP130; NOPP140; NS5ATP13; P130 ;

Immunogen:

Human P130.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A