

# SYT1 Antibody



PACO13553

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:****Protein Background:****Gene ID:**

SYT1

**Uniprot**

P21579

**Synonyms:**

v-rel reticuloendotheliosis viral oncogene homolog A (avian); RELA; MGC131774; NFKB3; p65 ;

**Immunogen:**

Human p65.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A