

PABPC1 Antibody



PACO13555

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

PABPC1

Uniprot

P11940

Synonyms:

poly(A) binding protein, cytoplasmic 1; PABPC1; PAB1; PABP; PABP1; PABPC2; PABPL1 ;

Immunogen:

Human PABPC1, PABP.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A