

# PACS2 Antibody



PACO13556

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PACS2

**Uniprot**

Q86VP3

**Synonyms:**

phosphofurin acid, c cluster sorting protein 2; PACS2; FLJ25488; KIAA0602; PACS1L ;

**Immunogen:**

Human PACS2-Specific.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A