

# CPB1 Antibody



PACO13559

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CPB1

**Uniprot**

P15086

**Synonyms:**

carboxypeptidase B1 (tissue); CPB1; DKFZp779K1333; PASP; PCPB ;

**Immunogen:**

Human PASP.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A