

PLA2G12B Antibody



PACO13566

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

PLA2G12B

Uniprot

Q9BX93

Synonyms:

22kd;

Immunogen:

Human PLA2G12B-Specific.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.320° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A