

POU3F1 Antibody



PACO13569

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

POU3F1

Uniprot

Q03052

Synonyms:

POU class 3 homeobox 1; POU3F1; OCT6; OTF6; SCIP ;

Immunogen:

Human POU3F1-Specific.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A