

PPP1R11 Antibody



PACO13572

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PPP1R11

Uniprot

O60927

Synonyms:

protein phosphatase 1, regulatory (inhibitor) subunit 11; PPP1R11; HCG-V; HCGV; MGC125741; MGC125742; MGC125743; TCTE5; TCTEX5 ;

Immunogen:

Human PPP1R11-Specific.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A