

IL18 Antibody



PACO13574

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

IL18

Uniprot

Q14116

Synonyms:

interleukin 18 (interferon-gamma-inducing factor); IL18; IGIF; IL-18; IL-1g; IL1F4; MGC12320 ;

Immunogen:

Human Pre-IL18.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A