

# RBM8A Antibody



PACO13584

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

RBM8A

**Uniprot**

Q9Y5S9

**Synonyms:**

RNA binding motif protein 8A; RBM8A; BOV-1A; BOV-1B; BOV-1C; MDS014; RBM8; RBM8B; Y14; ZNRP; ZRNP1 ;

**Immunogen:**

Human RBM8A, Y14.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A