

PACO13586

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

RNLS

Uniprot

Q5VYX0

Synonyms:

chromosome 10 open reading frame 59; C10orf59; FLJ11218; RENALASE ; C10orf59; RNLS; MAO-C

Immunogen:

Human RENALASE.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A