

PACO13590

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

RGN

Uniprot

Q15493

Synonyms:

regucalcin (senescence marker protein-30); RGN; RC; SMP30 ;

Immunogen:

Human RGN, SMP30.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A