

RIMS4 Antibody



PACO13593

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

RIMS4

Uniprot

Q9H426

Synonyms:

regulating synaptic membrane exocytosis 4; RIMS4; C20orf190; RIM4; dJ781B1.3 ;

Immunogen:

Human RIMS4-Specific.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A