

RREB1 Antibody



PACO13597

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

RREB1

Uniprot

Q92766

Synonyms:

ras responsive element binding protein 1; RREB1; FINB; HNT; LZ321; RREB-1; Zep-1 ;

Immunogen:

Human RREB1-Specific.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A