

TCEA1 Antibody



PACO13613

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TCEA1

Uniprot

P23193

Synonyms:

transcription elongation factor A (SII), 1; TCEA1; GTF2S; SII; TCEA; TF2S; TFIS ;

Immunogen:

Human SII-TFIS, TCEA1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A