

PACO13618

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

GINS4

Uniprot

Q9BRT9

Synonyms:

GINS complex subunit 4 (Sld5 homolog); GINS4; MGC14799; SLD5 ;

Immunogen:

Human SLD5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A