

SLIT2 Antibody



PACO13619

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SLIT2

Uniprot

O94813

Synonyms:

slit homolog 2 (Drosophila); SLIT2; FLJ14420; SLIL3; Slit-2 ;

Immunogen:

Human SLIT2-Specific.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A