

GUSBP1 Antibody



PACO13620

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

GUSBP1

Uniprot

Q15486

Synonyms:

glucuronidase, beta pseudogene; SMA4; FLJ36702; FLJ39772; MGC22265; MGC60382; SMA3; b55C20.2 ;

Immunogen:

Human SMA4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A