

SPTBN1 Antibody



PACO13625

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

SPTBN1

Uniprot

Q01082

Synonyms:

spectrin, beta, non-erythrocytic 1; SPTBN1; ELF; SPTB2; betaSp11 ;

Immunogen:

Human SPTBN1-Specific.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A