

TGOLN2 Antibody



PACO13636

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TGOLN2

Uniprot

O43493

Synonyms:

trans-golgi network protein 2; TGOLN2; MGC14722; TGN38; TGN46; TGN48; TGN51; TTN2 ;

Immunogen:

Human TGOLN2, TGN46.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A