

TOMM70A Antibody



PACO13647

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

TOMM70

Uniprot

O94826

Synonyms:

translocase of outer mitochondrial membrane 70 homolog A (S. cerevisiae); TOMM70A; FLJ90470 ; KIAA0719, TOM70

Immunogen:

Human TOM70.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A