

# THBS1 Antibody



PACO13651

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

THBS1

**Uniprot**

P07996

**Synonyms:**

thrombospondin 1; THBS1; THBS; THBS-1; TSP; TSP-1; TSP1 ;

**Immunogen:**

Human THBS1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A