

NXNL1 Antibody



PACO13654

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NXNL1

Uniprot

Q96CM4

Synonyms:

nucleoredoxin-like 1; NXNL1; RDCVF; TXNL6 ;

Immunogen:

Human TXNL6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A