

UBASH3A Antibody



PACO13656

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

UBASH3A

Uniprot

P57075

Synonyms:

ubiquitin associated and SH3 domain containing, A; UBASH3A; CLIP4; STS-2; TULA ;

Immunogen:

Human UBASH3A-Specific.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A