

UNC5D Antibody



PACO13662

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

UNC5D

Uniprot

Q6UXZ4

Synonyms:

unc-5 homolog D (C. elegans); UNC5D; FLJ16019; KIAA1777; PRO34692; Unc5h4 ;

Immunogen:

Human UNC5D-Specific.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A