

# ZNRF2 Antibody



PACO13669

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:****Gene ID:**

ZNRF2

**Uniprot**

Q8NHG8

**Synonyms:**

24kd;

**Immunogen:**

Human ZNRF2-Specific.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.320&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A