

RBP1 Antibody



PACO13671

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

RBP1

Uniprot

P09455

Synonyms:

CRABP-I; CRBP; CRBP1; CRBPI; RBPC

Immunogen:

Human RBP1 fusion protein.

Storage:

PBS with 0.02% sodium azide and 50% glycerol pH 7.320°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A