

CUL7 Antibody



PACO13676

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CUL7

Uniprot

Q14999

Synonyms:

cullin 7; CUL7; KIAA0076; dJ20C7.5 ;

Immunogen:

Human CUL7.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A