

DLK2 Antibody



PACO13677

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

DLK2

Uniprot

Q6UY11

Synonyms:

delta-like 2 homolog (Drosophila); DLK2; EGFL9; MGC111055; MGC2487 ;

Immunogen:

Human DLK2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A