

# HORMAD1 Antibody



PACO13690

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

HORMAD1

**Uniprot**

Q86X24

**Synonyms:**

HORMA domain containing 1; HORMAD1; CT46; DKFZp434A1315; NOHMA; RP11-363I22.1 ;

**Immunogen:**

Human HORMAD1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A