

NCALD Antibody



PACO13694

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NCALD

Uniprot

P61601

Synonyms:

neurocalcin delta; NCALD; MGC33870; MGC74858 ;

Immunogen:

Human NCALD.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A