

OLFM2 Antibody



PACO13705

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

OLFM2

Uniprot

O95897

Synonyms:

olfactomedin 2; OLFM2; NOE2; NOELIN2; NOELIN2_V1; OlfC ;

Immunogen:

Human OLFM2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A