

OLFM3 Antibody



PACO13706

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

OLFM3

Uniprot

Q96PB7

Synonyms:

olfactomedin 3; OLFM3; NOE3; NOELIN3; NOELIN3_V1; NOELIN3_V2; NOELIN3_V3; NOELIN3_V4; NOELIN3_V5; NOELIN3_V6; OPTIMEDIN ;

Immunogen:

Human OLFM3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A